



## Product Information

Dulbecco's Modified Eagle Medium (DMEM) High Glucose (4.5 g/L), w/o L-Glutamine

Catalog no: GBDM05 / 05F

### Specification

Appearance	Red colored, clear solution
pH	7.00 -7.60
Storage & Shelf Life	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping Conditions	Ambient

### Formulation

Components	Concentration mg/L	Components	Concentration mg/L
<b>Vitamins:</b>		<b>Amino Acids:</b>	
Choline chloride	4.00	Glycine	30.00
D-Calcium Pantothenate	4.00	L-Arginine HCl	84.00
Folic Acid	4.00	L-Cystine 2 HCl	62.57
myo-Inositol	7.20	L-Histidine HCl H <sub>2</sub> O	42.00
Nicotinamide	4.00	L-Isoleucine	105.00
Pyridoxal HCl	4.00	L-Leucine	105.00
Riboflavin	0.40	L-Lysine HCl	146.00
Thiamine HCl	4.00	L-Methionine	30.00
<b>Inorganic Salts:</b>		L-Phenylalanine	66.00
CaCl <sub>2</sub> 2 H <sub>2</sub> O	265.00	L-Serine	42.00
Fe(NO <sub>3</sub> ) <sub>3</sub> 9 H <sub>2</sub> O	0.10	L-Threonine	95.00
KCl	400.00	L-Tryptophan	16.00
MgSO <sub>4</sub> 7H <sub>2</sub> O	200.00	L-Tyrosine 2 Na 2H <sub>2</sub> O	103.79
NaCl	6400.00	L-Valine	94.00
NaHCO <sub>3</sub>	3700.00		
NaH <sub>2</sub> PO <sub>4</sub> 2H <sub>2</sub> O	141.73		
<b>Other Components:</b>			
D-Glucose	4500.00		
Phenol Red Sodium Salt	15.90		

### Need Help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genixisbiotech.com](mailto:info@genixisbiotech.com)

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.