



## Product Information

Leibovitz's L-15 Medium, w/o L-Glutamine

Catalog no: GBLL01 / 01F

### Specification

Appearance	Red Colored, Clear solution
pH	8.00 -8.60
Storage & Shelf Life	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping Conditions	Ambient

### Formulation

Components	Concentration mg/l	Components	Concentration mg/l
<b>Vitamins:</b>		<b>Amino Acids:</b>	
Choline Chloride	1.00	Glycine	200.00
D-Calcium Pantothenate	1.00	L-Alanine	225.00
Folic Acid	1.00	L-Arginine	500.00
myo-Inositol	2.00	L-Asparagine	250.00
Nicotinamide	1.00	L-Cysteine	120.00
Pyridoxine HCl	1.00	L-Histidine	250.00
Riboflavin	0.10	L-Isoleucine	125.00
Thiamine Monophosphate HCl	1.00	L-Leucine	125.00
<b>Inorganic Salts:</b>		L-Lysine HCl	93.70
CaCl <sub>2</sub> 2 H <sub>2</sub> O	185.00	L-Methionine	75.00
KCl	400.00	L-Phenylalanine	125.00
KH <sub>2</sub> PO <sub>4</sub> (anhydrous)	60.00	L-Serine	200.00
MgCl 6 H <sub>2</sub> O	200.00	L-Threonine	300.00
MgSO <sub>4</sub> (anhydrous)	97.67	L-Tryptophan	20.00
NaCl	8000.00	L-Tyrosine	300.00
Na <sub>2</sub> HPO <sub>4</sub> (anhydrous)	190.00	L-Valine	100.00
<b>Other Components:</b>			
D-Galactose	900.00		
Phenol Red Sodium Salt	11.00		
Sodium Pyruvate	550.00		

### Need Help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genexisbiotech.com](mailto:info@genexisbiotech.com)

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.