



Product Information

MEM Vitamin Solution (100X)

Catalog no: GBMVS01/01F

General Information

MEM Vitamin Solution is a growth supplement that is added to cell culture medium to support cell viability and development. The Minimum Essential Medium (MEM) Vitamin Solution is made to have 100X of vitamins (MEM).

Specification

Appearance	Clear pale yellow liquid
Storage and shelf life	Store at +2°C to +25°C upon arrival. Once opened store at +2°C to +25°C and use within 24 months.
Shipping conditions	Ambient

Formulation

Components	Concentration (mg/L)
Vitamins	Vitamins
Choline chloride	100.0
D-Calcium pantothenate	100.0
Folic Acid	100.0
Nicotinamide	100.0
Pyridoxal hydrochloride	100.0
Riboflavin	10.0
Thiamine hydrochloride	100.0
i-Inositol	200.0
Inorganic Salts	Inorganic Salts
Sodium Chloride (NaCl)	8500.0

Need Help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.