



## Product Information

Nuclease-Free Water (Non DEPC Treated)

Catalog no: GBNFWND02/02F

### General Information

Nuclease-free Deionized water that has not undergone DEPC treatment has been filtered into the final bottle and autoclaved. It doesn't need to be prepared, mixed, or autoclaved before use. The endonuclease, exonuclease, and RNase activity of this nuclease-free water is thoroughly tested for contamination.

### Specification

Appearance	Clear liquid
Storage and shelf life	Store at 4°C.
Shipping conditions	Ambient

### Need Help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genexisbiotech.com](mailto:info@genexisbiotech.com)

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.