



## Product Information

RPMI 1640, w/o L-Glutamine, w/o Phenol Red

Catalog no: GBRP04 / 04F

### Specification

Appearance	Clear solution
pH	7.00 -7.60
Storage & Shelf Life	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping Conditions	Ambient

### Formulation

Components	Concentration mg/L	Components	Concentration mg/L
<b>Vitamins:</b>		<b>Amino Acids:</b>	
p-Amino Benzoic Acid	1.00	Glycine	10.00
D-Biotin	0.20	L-Arginine	200.00
Choline Chloride	3.00	L-Asparagine H <sub>2</sub> O	56.82
D-Calcium	0.25	L-Aspartic Acid	20.00
Pantothenate		L-Cystine 2HCl	65.20
Folic Acid	1.00	L-Glutamic Acid	20.00
myo-Inositol	35.00	L-Histidine HCl H <sub>2</sub> O	20.27
Nicotinamide	1.00	L-Hydroxy-L-Proline	20.00
Pyridoxine HCl	1.00	L-Isoleucine	50.00
Riboflavin	0.20	L-Leucine	50.00
Thiamine HCl	1.00	L-Lysine HCl	40.00
Vitamin B12	0.005	L-Methionine	15.00
<b>Inorganic Salts:</b>		L-Phenylalanine	15.00
Ca(NO <sub>3</sub> ) <sub>2</sub> 4 H <sub>2</sub> O	100.00	L-Proline	20.00
KCl	400.00	L-Serine	30.00
MgSO <sub>4</sub> 7H <sub>2</sub> O	100.00	L-Threonine	20.00
NaCl	6000.00	L-Tryptophan	5.00
NaHCO <sub>3</sub>	2000.00	L-Tyrosine 2Na 2H <sub>2</sub> O	28.83
Na <sub>2</sub> HPO <sub>4</sub>	800.00	L-Valine	20.00
<b>Other Components:</b>			
D-Glucose	2000.00		
L-Glutathione Reduced	1.00		

### Need Help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genexisbiotech.com](mailto:info@genexisbiotech.com)

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.