



Product Information

Water, suitable for cell culture, sterile-filtered

Catalog no: GBW01/01F

General Information

High-quality water for use in cell culture applications that has been distilled, deionized, and disinfected using UV light. It may be used for a number of tasks, including preparing cell culture media, buffers, and laboratory chemicals, washing plastic, glass, and glassware, and cleaning bioreactors. It can also be used for techniques like chromatography (HPLC) and spectrometry.

Specification

Appearance	Clear colorless liquid
Storage and shelf life	Store at +2°C to +25°C upon arrival. Once opened store at +2°C to +25°C and use within 24 months.
Shipping conditions	Ambient

Need Help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.