

Product Information

Tris-Glycine SDS Gel Running Buffer (5X)

Catalogue Number: GBPTGSB02

General Information

The Tris-Glycine SDS Gel Running Buffer (5X) is formulated for SDS-PAGE applications, providing a stable pH and ionic conditions essential for protein electrophoresis under denaturing conditions. The 5X concentrated format allows for convenient dilution to a 1X working solution, ensuring consistent results in laboratory workflows.

Product Specification

Appearance	: Clear, colorless solution
Storage & Shelf Life	: Store at Room temperature and use within 12 months.
Shipping Conditions	: Ambient
pH	: 8.3
Concentration	: 5X

Formulation

Components	Concentration g/ml
Tris base	7.571
Glycine	36.1
SDS (Sodium Dodecyl Sulfate)	2.5

For Research Use Only

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com

Disclaimer: User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.