

## Product Information

### Ham's F-10, with L-Glutamine (Ham's F-10 Nutrient Mix)

Catalogue Number: GBHF03/03F

## Product Specification

Appearance : Clear red orange solution

pH : 7.00 -7.60

Storage & Shelf Life : Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.

Shipping Conditions: Ambient

## Formulation

Components	Concentration mg/L
<b>Vitamins:</b>	
D-Biotin	0.024
Choline chloride	0.698
D-Calcium Pantothenate	0.715
Folic Acid	1.32
myo-Inositol	0.541
Nicotinamide	0.615
Pyridoxine HCl	0.206
Riboflavin	0.376
Thiamine HCl	1.00
Vitamin B12	1.36
<b>Inorganic Salts:</b>	
CaCl <sub>2</sub> 2 H <sub>2</sub> O	74.10
CuSO <sub>4</sub> 5 H <sub>2</sub> O	0.0025
FeSO <sub>4</sub> 7 H <sub>2</sub> O	0.834
MgSO <sub>4</sub>	74.64
KCl	285.00
KH <sub>2</sub> PO <sub>4</sub>	83.00
NaHCO <sub>3</sub>	1200.00
NaCl	7400.00
Na <sub>2</sub> HPO <sub>4</sub>	153.70
ZnSO <sub>4</sub> 7 H <sub>2</sub> O	0.0288
<b>Other Components:</b>	
D-Glucose	1100.00
Hypoxanthine	4.08
Lipoic Acid	0.21
Phenol Red Sodium Salt	1.3
Sodium Pyruvate	110.00
Thymidine	0.73

Components	Concentration mg/L
<b>Amino Acids:</b>	
L-Alanine	8.91
L-Arginine HCl	211.00
L-Asparagine H <sub>2</sub> O	15.01
L-Aspartic Acid	13.30
L-Cysteine HCl H <sub>2</sub> O	35.13
Glycine	7.51
L-Glutamine	146.00
L-Glutamic Acid	14.70
L-Histidine HCl H <sub>2</sub> O	21.00
L-Isoleucine	2.60
L-Leucine	13.10
L-Lysine HCl	29.30
L-Methionine	4.48
L-Phenylalanine	4.96
L-Proline	11.50
L-Serine	10.50
L-Threonine	3.57
L-Tryptophan	0.60
L-Tyrosine 2 Na 2H <sub>2</sub> O	2.61
L-Valine	3.50

This product is for research use only.

## Need help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genexisbiotech.com](mailto:info@genexisbiotech.com)

Disclaimer: User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.