

Product Information

Leibovitz's L-15 Medium, w/o L-Glutamine

Catalogue Number: GBLL01 / 01F

Product Specification

Appearance : Red Colored, Clear solution

pH : 8.00 -8.60

Storage & Shelf Life : Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.

Shipping Conditions: Ambient

Formulation

Components	Concentration mg/L
Vitamins:	
Choline chloride	1.00
D-Calcium Pantothenate	1.00
Folic Acid	1.00
myo-Inositol	2.00
Nicotinamide	1.00
Pyridoxal HCl	1.00
Riboflavin	0.10
Thiamine Monophosphate HCl	1.00
Inorganic Salts:	
CaCl ₂ 2H ₂ O	185.00
KCl	400.00
KH ₂ PO ₄ (anhydrous)	60.00
MgCl 6 H ₂ O	200.00
MgSO ₄ (anhydrous)	97.67
NaCl	8000.00
Na ₂ HPO ₄ (anhydrous)	190.00
Other Components:	
D-Galactose	900.00
Phenol Red Sodium Salt	11.00
Sodium Pyruvate	550.00

Components	Concentration mg/L
Amino Acids:	
Glycine	200.00
L-Alanine	225.00
L-Arginine	500.00
L-Asparagine	250.00
L-Cystine	120.00
L-Histidine	250.00
L-Isoleucine	125.00
L-Leucine	125.00
L-Lysine HCl	93.70
L-Methionine	75.00
L-Phenylalanine	125.00
L-Serine	200.00
L-Threonine	300.00
L-Tryptophan	20.00
L-Tyrosine	300.00
L-Valine	100.00

This product is for research use only.

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com