

## Product Information

### McCoy's 5A Media with L-Glutamine

Catalogue Number: GBMCM01

## Product Specification

Appearance : Clear red orange solution

pH : 7.00 -7.60

Storage & Shelf Life : Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.

Shipping Conditions: Ambient

## Formulation

Components	Concentration mg/L
<b>Vitamins:</b>	
4-Aminobenzoic Acid	1.00
Ascorbic Acid	0.56
D-Biotin	0.20
Choline chloride	5.00
D-Calcium Pantothenate	0.20
Folic Acid	10.00
myo-Inositol	36.00
Nicotinamide	0.50
Nicotinic Acid	0.50
Pyridoxal HCl	0.50
Pyridoxine HCl	0.50
Riboflavin	0.20
Thiamine HCl	0.20
Vitamin B12	2.00
<b>Inorganic Salts:</b>	
CaCl <sub>2</sub> 2 H <sub>2</sub> O	132.00
MgSO <sub>4</sub>	97.67
KCl	400.00
NaHCO <sub>3</sub>	2200.00
NaCl	6460.00
Na <sub>2</sub> HPO <sub>4</sub> (anhydrous)	504.00
<b>Other Components:</b>	
D-Glucose	3000.00
L-Glutathione Reduced	0.50
Meat Peptone	600.00
Phenol Red Sodium Salt	11.00

Components	Concentration mg/L
<b>Amino Acids:</b>	
L-Alanine	13.36
L-Arginine HCl	42.14
L-Asparagine H <sub>2</sub> O	45.03
L-Aspartic Acid	19.97
L-Cysteine HCl H <sub>2</sub> O	31.75
L-Glutamic Acid	22.07
L-Glutamine	219.15
Glycine	7.51
L-Histidine HCl H <sub>2</sub> O	20.96
L-Hydroxyproline	19.67
L-Isoleucine	39.36
L-Leucine	39.36
L-Lysine HCl	36.54
L-Methionine	14.92
L-Phenylalanine	16.52
L-Proline	17.27
L-Serine	26.28
L-Threonine	17.87
L-Tryptophan	3.06
L-Tyrosine 2 Na	26.10
L-Valine	17.57

This product is for research use only.

## Need help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genexisbiotech.com](mailto:info@genexisbiotech.com)

Disclaimer: User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.