

Product Information

PBS (10X) Endotoxin Free

Catalogue Number: GBPPB01/01F

General Information

10X PBS, Endotoxin-Free is a high-quality, sterile solution designed for use in sensitive cell culture, molecular biology, and immunology experiments. PBS is isotonic and non-toxic to cells, making it ideal for maintaining the pH and osmotic balance in various biological applications. The 10X concentrate can be diluted to 1X with endotoxin-free water for routine use, ensuring convenience and consistent results in endotoxin-sensitive assays.

Product Specification

Appearance : Clear, colorless solution

pH : 7.4

Storage & Shelf Life : Store at Room temperature and use within 12 months

Endotoxin activity : <0.005 EU/mL

Shipping Conditions : Ambient

Formulation

Components	Concentration g/L
NaCl	80.1
KCl	2
Na ₂ HPO ₄	14.4
KH ₂ PO ₄	2.4

This product is for research use only.

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com

Disclaimer: User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.